



An End-To-End Precision Medicine Approach for Matching Cancer Patients to Clinical Trials

Alexei Protopopov, Jia Xu, Claudia S. Huettner, Meaghan Russell, Scott Schell, Stephen Lyle

KEW Group, Inc., Cambridge, MA; Department of Cancer Biology, University of Massachusetts Medical School, Worcester, MA

INTRODUCTION

Precise matching of a patient's tumor genomic profile with appropriate targeted agents is the goal of precision medicine, and the salutary benefits of these targeted therapies upon survival has been well documented in the literature. Using a 413+ cancer-related actionable gene, NGS sequencing panel, KEW CancerPlex®, we analyzed a variety of solid tumors to identify mutations associated with predicted response to approved and investigational targeted therapy agents.

METHODS

Tumor profiling was conducted in a CAP/CLIA-certified laboratory, licensed by MA and 48 other states (KEW Group Inc.). DNA was extracted from FFPE tissue sections, slides, cell blocks, from FNAs or effusions, or cell pellets, followed by hybrid-capture next-gen sequencing. Rapid sequencing runs were employed to generate at least 200x depth, and mutational analysis was performed using the KEW Clinical Genomics Analytical pipeline. The coding regions and portions of the introns of 413+ genes were sequenced. The assay simultaneously surveyed multiple classes of genomic abnormalities including single nucleotide substitutions (SNP), small insertions/deletions (indels), copy number alterations (CNV), and translocations.

Variant characterization and annotation were conducted for a MAF cut-off of 10%. Genotype-based personalized molecular modeling was performed to characterize tumors, and response or resistance to FDA-approved drugs was predicted. Matches were identified by screening 260,000 clinical trials, including approximately sixty thousands active or recruiting studies.

TECHNICAL PERFORMANCE - CancerPlex® v3

Analytic sensitivity for SNP calls is 98.9%-100% at 95% confidence interval, specificity is 99.99%; for indels sensitivity is 99.8%-100% at 95% confidence interval and specificity is 99.99%. Analytic sensitivity for translocations of the ALK, RET, and ROS1 genes is 100% for tumor purity as low as 10%.



Testing and reporting turnaround time of 7-10 days.

CancerPlex Test requires only a small portion of the biopsy (40 micron section or equivalency), and it can technically be performed with a very low amount of input DNA. The CancerPlex requirement of tissue is determined by biological considerations to obtain a representative sampling of the tumor.

RESULTS

200 cases, representing a variety of tumor types were analyzed using CancerPlex® v2 Assay. Tumors averaged 6.65 somatic variants (range: 1-18), of which 2.48 variants (range: 0-9) were deemed clinically actionable (93%). Tumor response, or resistance markers were identified for FDA-approved therapies (mean: 1.2; range: 0-5) and clinical trials (mean: 7.98 range: 0-24).


Additional 50 randomly selected clinical cases were characterized using the upgraded CancerPlex® v3 Assay that provides about 2-fold deeper sequence coverage (500x) for the same gene panel. Overall fifty seven alterations were found actionable, those alterations are assigned to 47 cases (94%). None of the samples failed to yield a report.



ACTIONABILITY

Our study shows that the actionability rate is not significantly affected by an increase of sequencing coverage. The major reason is the arbitrary cut-off of 10% MAF. The actionability rate could be improved if one were to consider including clonal alterations assigned to only a subset of the entire tumor cell population.

Another factor impacting actionability may be the gene panel size. Published validation studies show that hot-spot analysis (small panels of ~100 cancer genes) misses up to 20% of clinically important mutations in some tumors (Drilon et al., 2015). To test this further, we took a subset of CancerPlex genes to match the content of a previously published smaller panel of 315 genes (Frampton et al., 2015), and compared the actionability. The larger panel provides critically important additions (see Figures).



CONCLUSION

Larger NGS panel size increases the likelihood of identifying clinically actionable findings, including FDA-approved therapies missed by NGS panels with fewer genes.

CancerPlex Current Gene List: www.kewgroup.com/sites/default/files/uploads/downloads/kewcancerplexfactsheet.pdf

ABL1	ABL2 (ARG, BACH)	ACVR1B	ACVR2A	AFF1	AKT1	AKT2	AKT3	ALK	AMER1 (FAM123B)	APC	AR	ARAF	ARFRP1 (ARP1)	ARID1A	ARID2	ARNT (HIF-1β)	ASCL4	ASXL1	ATM	ATR	ATRX	AURKA	AURKB	AURKC	AXIN2	AXL (UFO)	BAP1	BARD1	BCL10																																																																														
BCL11A (ZNF856)	BCL11B	BCL2	BCL2L1 (BCL-X)	BCL2L2 (BCL-X)	BCL6	BCL9	BCOR	BCORL1	BCORL2 (IP2)	BIRC3 (Aip1)	BLM	BLNK	BMPR1A	BRAF	BRCA1	BRCA2	BRIP1 (FANCI)	BTBK	BUB1B	C11orf30 (EMSY)	CARD11	CASP8	CBFB	CBL	CCND1 (Cyclin D1)	CCND2 (Cyclin-D2)	CCND3 (Cyclin D3)	CCNE1 (Cyclin E1)	CD70 (TNFSF7)																																																																														
CD79A	CD79B	CD373 (HRPT1)	CDH1 (E-Cadherin)	CDH11 (Catherin 11)	CDH2 (N-Cadherin)	CDK12	CDK4	CDKN1A	CDKN2A (ARF, INK4)	CDKN2B (INK4B)	CDKN2C (INK4C)	CEBPA	CHEK1	CIC	CREBBP	CRKL	CRF2	CSF1R	CSMD3	CTCF	CTNNA1 (catenin)	CTNNA2 (catenin)	CTNNA3 (catenin)	CTNNA4 (catenin)	CTNNA5 (catenin)	CTNNA6 (catenin)	CTNNA7 (catenin)	CTNNA8 (catenin)	CTNNA9 (catenin)																																																																														
CYLD	DAXX	DCC	DDIT3 (CHOP)	DDIT3 (CHOP)	DDR1	DDR2	DICER1	DNMT3A	DNMT3B	DNMT3C	DST	EDNRB	EGFR	EMIL4	EP300 (P300)	EPHA3	EPHA5	EPHA7	EPHB1	EPHA4	EPHB6	EPHB7	EPHB8 (HER2)	EPHB9 (HER3)	EPHB10 (HER4)	ERCC1 (RAD10)	ERCC2 (XPD)	ERCC3 (XPB, RAD25)	ERCC4 (XPF)	ERCC5																																																																													
ESR1	ETS1	ETV1	ETV4	EXT1	EXT2	FOXO1	FOXO3	FANCA	FANCC	FANCD1	FANCF	FANCG	FANCL	FAS (APO-1)	FAT1	FBXW7 (BAF250)	FGF10	FGF14	FGF19	FGF23	FGF3	FGF4	FGF5	FGF6	FGF7	FGF8	FGF9	FGFR1	FGFR2	FGFR3	FGFR4	FLCN																																																																											
FLT1	FLT1 (VEGFR1)	FLT3 (STK1)	FLT4 (VEGFR-3)	FN1	FOXO1	FOXO3	FOXO4	GATA1	GATA2	GATA3	GATA4	GATA5	GATA6	GATA7	GATA8	GATA9	GATA10	GATA11	GATA12	GATA13	GATA14	GATA15	GATA16	GATA17	GATA18	GATA19	GATA20	GATA21	GATA22	GATA23	GATA24	GATA25	GATA26	GATA27	GATA28	GATA29	GATA30	GATA31	GATA32	GATA33	GATA34	GATA35	GATA36	GATA37	GATA38	GATA39	GATA40	GATA41	GATA42	GATA43	GATA44	GATA45	GATA46	GATA47	GATA48	GATA49	GATA50	GATA51	GATA52	GATA53	GATA54	GATA55	GATA56	GATA57	GATA58	GATA59	GATA60	GATA61	GATA62	GATA63	GATA64	GATA65	GATA66	GATA67	GATA68	GATA69	GATA70	GATA71	GATA72	GATA73	GATA74	GATA75	GATA76	GATA77	GATA78	GATA79	GATA80	GATA81	GATA82	GATA83	GATA84	GATA85	GATA86	GATA87	GATA88	GATA89	GATA90	GATA91	GATA92	GATA93	GATA94	GATA95	GATA96	GATA97	GATA98	GATA99	GATA100